



## Model SM161G Glass Door Pharmacy Refrigerator



### Shoreline SM161G Glass Door Pharmacy Refrigerator

- External : 835mmH x 600mmW x 615mmD
- Internal : 667mmH x 441mmD x 439mmD (*including rear step*)
- 141 litres Capacity
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display with ShoreTouch™ controls
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Self-Closing Door with Lock and 2 keys
- High Efficiency - Thick Wall Insulation
- Adjustable Shelving x 3
- Interior Light
- CFC & Ammonia free
- RPSGB, NHS, WHO, RCVS compliant
- 220-230v 50Hz
- Max. Ambient +38°C
- Made in England
- 2 Years Parts & Labour Warranty (UK only)
- Supplied with Certificate of Conformity

(December 2010)