



Model SM161 Solid Door Pharmacy Refrigerator



Shoreline SM161 Solid Door Pharmacy Refrigerator

- External : 835mmH x 600mmW x 615mmD
- Internal : 667mmH x 441mmD x 439mmD (including rear step)
- 141 litres Capacity
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display with ShoreTouch™ controls
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Self-Closing Door with Lock and 2 keys
- High Efficiency - Thick Wall Insulation
- Adjustable Shelving x 3
- CFC & Ammonia free
- RPSGB, NHS, WHO, RCVS compliant
- 220-230v 50Hz
- Max. Ambient +43°C
- Made in England
- 2 Years Parts & Labour Warranty (UK only)
- Supplied with Certificate of Conformity

(December 2010)

SHORELINE (UK) LTD

16-17 Martello Enterprise Centre, Courtwick Lane, Littlehampton, West Sussex BN17 7PA

Tel: 01903 733877 Fax: 01903 733891 www.shoreline-medical.co.uk ask@shoreline-medical.co.uk