



Model SM360G Glass Door Pharmacy Refrigerator



Shoreline SM360G Glass Door Pharmacy Refrigerator

1595mmH x 600mmW x 600mmD
352 litres Capacity
Temperature Range +2°C to +8°C
Suitable for Vaccine & Pharmaceutical Storage
Forced-Air Cooling for even temperature distribution
Digital LED Temperature Display
Audio / Visual Alarm on Temperature Deviation
MAX/MIN Temperature Memory
Automatic Defrosting
Door Lock and 2 keys
CFC & Ammonia free
Adjustable Shelving x 5
Interior Light
RPSGB, NHS, WHO, RCVS compliant
220-230v 50Hz
Max. Ambient +32°C
Made in England
2 Years Part & Labour Warranty (UK only)

(May 2010)

SHORELINE (UK) LTD

16-17 Martello Enterprise Centre, Courtwick Lane, Littlehampton, West Sussex BN17 7PA

Tel: 01903 733877 Fax: 01903 733891 www.shoreline-medical.co.uk ask@shoreline-medical.co.uk