



## Model SM541G Pharmacy Refrigerator



### Shoreline SM541G Glass Door Pharmacy Refrigerator

1745mmH x 750mmW x 700mmD

544 litres Capacity

Temperature Range +2°C to +8°C: Suitable for Vaccine & Pharmaceutical Storage

Forced-Air Cooling for even temperature distribution

Digital LED Temperature Display

Audio / Visual Alarm on Temperature Deviation

MAX/MIN Temperature Memory

Remote Alarm Terminals

Automatic Defrosting

Self-Closing Door with Lock and 2 keys

Castors (2 front locking)

High Efficiency - Thick Wall Insulation

1.4 kWhr/24hr

Adjustable Shelving x 5

Interior Light

CFC & Ammonia free

RPSGB, NHS, WHO, RCVS compliant

220-230v 50Hz

Max. Ambient +38°C

Made in England

2 Years Part & Labour Warranty (UK only)

(May 2010)

SHORELINE (UK) LTD

16-17 Martello Enterprise Centre, Courtwick Lane, Littlehampton, West Sussex BN17 7PA

Tel: 01903 733877 Fax: 01903 733891 [www.shoreline-medical.co.uk](http://www.shoreline-medical.co.uk) [ask@shoreline-medical.co.uk](mailto:ask@shoreline-medical.co.uk)