



## Model SM541 Pharmacy Refrigerator



### Shoreline SM541 Solid Door Pharmacy Refrigerator

1745mmH x 750mmW x 700mmD  
544 litres Capacity  
Temperature Range +2°C to +8°C: Suitable for Vaccine & Pharmaceutical Storage  
Forced-Air Cooling for even temperature distribution  
Digital LED Temperature Display  
Audio / Visual Alarm on Temperature Deviation  
MAX/MIN Temperature Memory  
Remote Alarm Terminals  
Automatic Defrosting  
Self-Closing Door with Lock and 2 keys  
Castors (2 front locking)  
High Efficiency - Thick Wall Insulation  
1.0 kWhr/24hr  
Adjustable Shelving x 5  
CFC & Ammonia free  
Meets RPSGB, NHS, WHO, RCVS performance criteria  
220-230v 50Hz  
Max. Ambient +43°C  
Made in England  
2 Years Part & Labour Warranty (UK only)

(May 2010)

SHORELINE (UK) LTD

16-17 Martello Enterprise Centre, Courtwick Lane, Littlehampton, West Sussex BN17 7PA

Tel: 01903 733877 Fax: 01903 733891 [www.shoreline-medical.co.uk](http://www.shoreline-medical.co.uk) [ask@shoreline-medical.co.uk](mailto:ask@shoreline-medical.co.uk)