



Model SM75 Solid Door Pharmacy Refrigerator



Shoreline SM75 Solid Door Pharmacy Refrigerator

- External Dimensions : 790mmH x 450mmW x 430mmD
- Internal Dimensions : 525mmH x 375mmW x 305mmD
- Capacity : Gross 75 litres / Net 72 litres
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display with ShoreTouch™ controls
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and 2 keys
- CFC & Ammonia free
- Adjustable Shelving x 2
- RPSGB, NHS, WHO, RCVS compliant
- 220-240v 50/60 Hz
- Max. Ambient +32°C
- Made in England
- 2 Years Parts & Labour Warranty (UK only)
- Supplied with Certificate of Conformity

(January 2011)

SHORELINE (UK) LTD

16-17 Martello Enterprise Centre, Courtwick Lane, Littlehampton, West Sussex BN17 7PA

Tel: 01903 733877 Fax: 01903 733891 www.shoreline-medical.co.uk ask@shoreline-medical.co.uk