



## ShoreLog™ Temperature Data Logger Kit



## ShoreLog™ Temperature Data Logger Kit

- Ensure your vaccine fridge temperatures are recorded 24/7 using the new ShoreLog Temperature Data Logger. Simply hang inside your fridge and download recorded temperatures in a report format once a week. Any deviations are detailed in against date & time interval via a simple print report menu.
- Kit comprises of 1 x Data Logger Capsule on coloured key ring tag, USB Communications Cable, Software CD for single PC, and User Guide Manual
- Ready for connection direct into your PC or Laptop computer
- Serial Port interface option (on request)
- Logging time intervals can be set between 1 minute and 255 minutes
- Extensive reporting options; Graphs, Lists, Deviations from set-point range
- Compatible with Windows XP or Vista, Windows 7.
- Battery Life: 10 years or 1 million readings maximum
- 12 months technical telephone support
- Additional Capsules available in Singles, or in Packs of 5, or 10

(August 2010)